



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

Retrograde Ejaculatory Urine/Semen

Test Code:

Specimen Type: Urine and Semen Samples

Tube Type / Collection Container: 2 Clean Containers Labeled A and B

Specimen Handling: Induce ejaculation by masturbation into container A, cleanse penis with water. Void urine into container B. Keep semen sample close to body temperature during transport. See Laboratory Services Customer Manual for complete details.

Specimen Storage: Body Temperature

Specimen Collection Notes: Deliver to lab within 1 hour of collection Monday - Friday 7:00 a.m. to 7:00 p.m. Saturday 7:00 a.m. to noon. Give patient a copy of collection instructions with container.

Department: Urinalysis

