



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

Antepartum Group B Streptococcus Screen by PCR (35-37 week screening)

Test Code: SSDL

Specimen/Specimen Collection: Vaginal/rectal specimen collected on Copan Venturi Transystem double swab. Specimens received on other swabs/transport devices will be rejected for testing.

Specimen Handling: If transported to the laboratory within 24 hours the specimen may be kept at room temperature. After 24 hours, the specimen must be refrigerated at 2-8°C. Maximum specimen age: 6 days at 2-8°C.

Department: Microbiology



CopanVenturi Transystem double swab