



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

ANGIOTENSIN 1 CONVERTING ENZYME

Test Code: ANGIO

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 2 ML

Specimen Handling: Centrifuge

Specimen Storage: Room Temperature

Specimen Collection Notes: Stop administration of Captopril, Enalapril, Lisinopril for 12 hrs. prior to collection.

Department: Reference Laboratory

Reference Lab Code: 010116

