

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

B. Pertussis Smear

Test Code: MCB

Specimen Type: Nasal Swab

Tube Type / Collection Container: 3-4 Slides from Nasal Swab

Specimen Volume: 3-4 slides

Specimen Handling: Air Dried Slides Do NOT FIX


Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: 3-4 Microscope slides of material from nasal

swab. Air dry. Do NOT Fix. Room Temp

Department: Microbiology / Reference Laboratory

Reference Lab Code: 086173



3-4 Air Dried Slides - DO NOT FIX