



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## **COPPER/SERUM or PLASMA**

**Test Code:** COPPE

**Specimen Type:** Serum or Plasma

**Tube Type / Collection Container:** Brick Color Top OR Dark Blue Trace Metal Tube ;EDTA

**Specimen Volume:** 1 ML

**Specimen Handling:** Centrifuge and pipette plasma or serum into plastic screw cap transport tube.

**Specimen Storage:** Ambient / Room Temp.

**Specimen Collection Notes:** Serum must be separated from cells w/in 45 mins of collection

**Department:** Reference Laboratory

**Reference Lab Code:** 001586



OR

