



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

B ANTHRACIS CULTURE

Test Code: ANTHR

Specimen Type: Whole Blood or Stool or Sputum or CSF

Tube Type / Collection Container: Swab or Blood Culture Tubes or Sterile Container

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Skin lesions submit swab with Stewarts media; Blood submit blood culture bottles; CSF, Stool and Sputum submit in sterile container. Transport ASAP. Specimens stable 24 hours at room temperature.

Department: Microbiology

