



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

FETAL FIBRONECTIN

Test Code: FETFB

Specimen Type: Swab

Tube Type / Collection Container: Rapid fFN Specimen Collection Kit (Hologic Inc)

Specimen Volume: Transport tube contains required media do not add liquid

Specimen Handling: Swab placed into the kit transport media

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Kit contains collection instructions

Department: Chemistry



