



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## GENTAMICIN, PEAK

**Test Code:** GENTP

**Specimen Type:** Serum

**Tube Type / Collection Container:** BRICK OR RED TOP TUBE

**Specimen Volume:** 2 ML Collect 60 to 90 minutes after dose

**Specimen Handling:** Centrifuge and pipette Serum into plastic screw cap transport tube within 1HR of collection.

**Specimen Storage:** Refrigerate 2-8 C Stability: 2 HRS at Room Temp and 7 Days Refrigerated.

**Department:** Chemistry



OR

