



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

GENTAMICIN,TROUGH

Test Code: GENTT

Specimen Type: Serum

Tube Type / Collection Container: BRICK OR RED TOP TUBE

Specimen Volume: 2 ML Collect immediately prior to dose and not longer than 30 minutes prior to dose.

Specimen Handling: Centrifuge and pipette Serum into plastic screw cap transport tube within 1 HR of collection.

Specimen Storage: Refrigerate 2-8 C Stability: 2 HRS at Room Temp and 7 Days Refrigerated.

Department: Chemistry



OR



