



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

GLUCOSE, 1 HR POST

(50 gram load)

Test Code: GL50

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 1 ML

Specimen Handling: Collect at ECH. Centrifuge. 12 HR Fast required. Collect Specimen 1 HR after oral dose is given. Patient must remain at the Hospital during the 1 HR testing period.

Specimen Storage: Refrigerate 2-8 C. Stability 48 HRS

Department: Chemistry



