

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

GLIADIN PROFILE –IgA,IgG,IgM Method: EIA

Test Code: MISC

Specimen Type: Serum

Tube Type / Collection Container: SST SERUM SEPARATOR

Specimen Volume: 1 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate

Specimen Collection Notes: Specimen must be free of Hemolysis, Lipemia

and Bacterial Contamination.

Department: Reference Laboratory

Reference Lab Code: 163402

