



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

PROSTATIC SPECIFIC ANTIGEN, TOTAL SCREEN

Test Code: PSASC

Specimen Type: serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 2 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Specimen stable 7 days refrigerated

Department: Chemistry

