



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

RAPID “A” STREP SCREEN, THROAT

Test Code: RSSA

Specimen Type: Serum

Tube Type / Collection Container: SWAB –STUART MEDIA (RAYON OR POLY)

Specimen Volume:

Specimen Handling: Break/ Crush the Ampule in bottom of swab.

Specimen Storage: Refrigerate 2-8 C Stable 4 HRS

Department: Hematology



Swab may vary in appearance slightly – Stuarts Media is Required

