



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

IGA IMMUNOGLOBULIN

Test Code: IGAQ

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 1 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: 12hr fasting required/serum stable 1 week refrigerated/min. 0.4ml pediatric

Department: Chemistry

