



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

MICROALBUMIN, URINE/RNDM

Test Code: MALBR

Specimen Type: Urine

Tube Type / Collection Container: Clean Container

Specimen Volume: 10 ML

Specimen Handling:

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Random or 24hr Specimen Accepted No preservative. Microalbumin/Creatinine Ratio reported with this test.

Department: Chemistry



Clean Container – Sterile Cup NOT Required