



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

ARSENIC Exposure Profile

Test Code: MISC

Specimen Type: Urine

Tube Type / Collection Container: 24 HR JUG OR Clean Container

Specimen Volume: Submit Entire 24HR Collection OR Random

Specimen Handling: 24 HR Urine or Randomly Voided

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Patient must avoid all seafood for 72hrs prior to collection / **Sampling time is the end of the work week for industrial exposure monitoring.** / Container: Plastic 24hr No preservative . Follow standard 24hr urine collection procedures. Must state 24hr volume.

Department: Reference Laboratory

Reference Lab Code: 007039



No Preservative OR

