

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

## CALCIUM, TIMED URINE

Test Code: CAU

Specimen Type: Urine

Tube Type / Collection Container: Clean Container

Specimen Volume: 10 ML

**Specimen Handling: Preservative: 10 ML 6N HCL** 

**Specimen Storage:** Ambient / Room Temp.

**Specimen Collection Notes: Only with physician approval**; Patient should be on a low calcium, neutral ash diet for 3 days prior to collection. Drugs affecting mineral metabolism should be withdrawn, if possible, 2-4 weeks prior to and during collection. / Container: Plastic w/10 ML 6N HCL;must state date and time of random collection. PH must be <4

**Department:** Reference Laboratory

**Reference Lab Code: 013706** 



Actual Cup may vary in Appearance - Sterile Cup is NOT Required