



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

CELL COUNT, CSF

Test Code: CELLC

Specimen Type: CSF

Tube Type / Collection Container: CSF - sterile plastic tube

Specimen Volume: 1 ML

Specimen Handling:

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Send to Lab within 1hr. Or Refrig.and send ASAP / Test includes WBC, RBC, Diff

Department: Hematology

