



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## CELL COUNT, JOINT FLUID

**Test Code:** CELJT

**Specimen Type:** Joint Fluid

**Tube Type / Collection Container:** Dark Green Sodium Heparin or  
Lavender K3 EDTA

**Specimen Volume:** 1 ML Specify Site-Source

**Specimen Handling:** Mix well

**Specimen Storage:** Refrigerate 2-8 C

**Specimen Collection Notes:** Send to Lab within 1hr. Or Refrig.and send  
ASAP / Test includes WBC, RBC, Diff

**Department:** Hematology

