



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

FOLATE RBC w/HCT

Test Code: FOLRC

Specimen Type: Whole Blood

Tube Type / Collection Container: 2 K3 EDTA

Specimen Volume: 2 Full K3EDTA tubes

Specimen Handling: 1 K3 EDTA Refrigerated 2-8 C, 1 K3 EDTA whole blood transferred to plastic screw top transport tube and Frozen.

Specimen Storage: Refrigerate 2-8 C K3 EDTA tube / Frozen Whole Blood

Department: Reference Laboratory

Reference Lab Code: 266015

