



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

Fatty acids, Free

Test Code: MISC

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 1 ML

Specimen Handling: Centrifuge

Specimen Storage: Frozen

Specimen Collection Notes: Sample should be collected in the early morning after a 12-hour fast. Clot, Separate and Freeze serum immediately.

Department: Reference Laboratory

Reference Lab Code: 081893