



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

JOINT FLUID ANALYSIS

Test Code: JTFLD

Specimen Type: Fluid

Tube Type / Collection Container: Dark Green; Sodium Heparin

Specimen Volume: 3 ML Specify Site - Source

Specimen Handling:

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Send to Lab within 1hr. Or Refrig. and send ASAP. Test includes WBC, RBC, Diff, Crystals, Glucose, Protein, Mucin Clot and RA

Department: Hematology



