



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

OCCULT BLOOD, 1-3 DETERMINATIONS

Test Code: HEMOC

Specimen Type: Feces

Tube Type / Collection Container: KIT (Hemoccult)

Specimen Volume: stool smear

Specimen Handling:


Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Follow Lab collection instructions and return to the lab as soon as possible, no longer than 14 days after collection. Do not get kit wet and maintain at room temperature before and after collection.

Department: Urinalysis



Front of Card



Back of Card