



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

SEMEN EXAM - POST VASECTOMY

Test Code: SEMPV

Specimen Type: Semen

Tube Type / Collection Container: Clean Container

Specimen Volume: n/a

Specimen Handling: Maintain specimen at body temperature by transporting close to body. E.g. inside a pocket.

Specimen Storage: Body Temperature

Specimen Collection Notes: Deliver to lab within 1 hour of collection
Monday - Friday 7:00 a.m. to 7:00 p.m. Saturday 7:00 a.m. to noon. Give patient a copy of collection instructions with container.

Department: Urinalysis



