



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

SODIUM, URINE 24 HOUR

Test Code: NA24U

Specimen Type: Urine

Tube Type / Collection Container: 24 hr Urine Container

Specimen Volume: Submit entire 24 HR urine collection

Specimen Handling: Submit entire 24 HR urine collection

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Refer to 24 hour urine collection procedure in this manual

Department: Chemistry



This Collection Requires NO Preservatives

