

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

VMA 24 Hour Urine -Vanillylmandelic Acid

Test Code: VMAU

Specimen Type: Urine

Tube Type / Collection Container: 24 hr Urine Container

Specimen Volume: 30 ML

Specimen Handling: Submit entire urine collection

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Record total volume; preservative - 30ML 6N HCl; See instructions for additional information / Specimen is stable for 5 days roomtmp w/no preserv.-12/23/04 Refer to 24 hour urine collection procedure in this manual

Department: Reference Laboratory

Reference Lab Code: 004143



30ML 6N HCl